



*Optimum quality  
at the fair price*



K1/2011

# *"MEDICA" Enterprise Price List for Dental Equipment*



## Dental unit

### ABC 04-012 (1 modification)

#### Chair-mounted



<b>ABC 04-012 (1modification) chair-mounted, with lower handpieces delivery system.</b> In set: block, patient chair with left plastic armrest, saliva ejector suction system for spittle equipment table for five equipments on pantograph arm, bottle water system, multifunctional pedal and table for instruments.	<b>4 130 €</b>
<b>ABC 04-012 (1 modification) chair-mounted, with upper handpieces delivery system</b>	<b>4 370 €</b>
<b>Autonomic drain system</b>	<b>140 €</b>
<b>Complekting operating light</b>	
Halogen dental light with two light intensity on pantograph arm	<b>223 €</b>
LED dental light DL60/DL66 with 1 LED, two light intensity on pantograph arm	<b>376 €</b>
LED dental light DL150/DL155 with 4 LED, three light intensity on pantograph arm	<b>400 €</b>
LED dental light with 6 LED, three light intensity on pantograph arm	<b>360 €</b>
<b>Variable equipment set for unit</b>	
Turbine connector	<b>133 €</b>
Turbine connector with light	<b>188 €</b>
Micromotor MV2.20 with connector	<b>160 €</b>
Air micromotor VORTICE (with turbo connection)	<b>300 €</b>
Micromotor brushless with connector DEFINITIVE SET	<b>350 €</b>
Micromotor brushless with LED, with connector DEFINITIVE LED SET	<b>435 €</b>
Technical micromotor 102 L with connector	<b>266 €</b>
Technical micromotor 108 E with connector	<b>226 €</b>
Multifunctional syringe – angular, 3 function	<b>70 €</b>
Multifunctional syringe – straight, 3 function	<b>115 €</b>
Multifunctional syringe – straight /angular, 6 function	<b>220 €</b>
Multifunctional syringe Minilight Fibra Optical – straight /angular, 7 function	<b>261 €</b>
Multifunctional syringe Minibright – straight /angular, 7 function	<b>232 €</b>
LED curing lamp	<b>480 €</b>
Ultrasonic scaler	<b>510 €</b>
Airpolishing for turbine coupling „Easyjet Pro“	<b>460 €</b>
Negotoscope	<b>130 €</b>
<b>Suction system with assistant table</b>	
Injector air suction system for dust	<b>270 €</b>
Vacuum suction system with automatic air/water separation “METASYS” (without vacuum pump)	<b>410 €</b>
<b>Additional dental equipment completing</b>	
Intra-oral camera with monitor and holder	<b>1 235 €</b>
Vacuum pump in box ABC 04-VS	<b>695 €</b>
Table for devices	<b>45 €</b>
Two metal armrests for patient chair (can be one armrest)	<b>185 €</b>
Digital control for micromotor DEFINITIVE	<b>376 €</b>

## Dental unit

### ABC 04-012 – Amber



<b>SIGER U 100 chair-mounted, with lower handpieces delivery system.</b> In set: block, patient chair with left plastic armrest, injector air suction system for spittle, equipment table for five equipments on pantograph arm, bottle water system, multifunctional pedal and table for instruments.	<b>4 830 €</b>
<b>SIGER U 100 chair-mounted, with upper handpieces delivery system</b>	<b>5 070 €</b>
<b>Autonomic drain system</b>	<b>140 €</b>
<b>Complekting operating light</b>	
Halogen dental light with two light intensity on pantograph arm	<b>223 €</b>
LED dental light DL60/DL66 with 1 LED, two light intensity on pantograph arm	<b>376 €</b>
LED dental light DL150/DL155 with 4 LED, three light intensity on pantograph arm	<b>400 €</b>
LED dental light with 6 LED, three light intensity on pantograph arm	<b>360 €</b>
<b>Variable equipment set for unit</b>	
Turbine connector	<b>133 €</b>
Turbine connector with light	<b>188 €</b>
Micromotor MV2.20 with connector	<b>160 €</b>
Air micromotor VORTICE (with turbo connection)	<b>300 €</b>
Micromotor brushless with connector DEFINITIVE SET	<b>350 €</b>
Micromotor brushless with LED, with connector DEFINITIVE LED SET	<b>435 €</b>
Technical micromotor 102 L with connector	<b>266 €</b>
Technical micromotor 108 E with connector	<b>226 €</b>
Multifunctional syringe – angular, 3 function	<b>70 €</b>
Multifunctional syringe – straight, 3 function	<b>115 €</b>
Multifunctional syringe – straight /angular, 6 function	<b>220 €</b>
Multifunctional syringe Minilight Fibra Optical – straight /angular, 7 function	<b>261 €</b>
Multifunctional syringe Minibright – straight /angular, 7 function	<b>232 €</b>
LED curing lamp	<b>480 €</b>
Ultrasonic scaler	<b>510 €</b>
Airpolishing for turbine coupling „Easyjet Pro“	<b>460 €</b>
Negotoscope	<b>130 €</b>
<b>Suction system with assistant table</b>	
Injector air suction system for dust	<b>270 €</b>
Vacuum suction system with automatic air/water separation “METASYS” (without vacuum pump)	<b>410 €</b>
<b>Additional dental equipment completing</b>	
Intra-oral camera with monitor and holder	<b>1 235 €</b>
Vacuum pump in box ABC 04-VS	<b>695 €</b>
Table for devices	<b>45 €</b>
Two metal armrests for patient chair (can be one armrest)	<b>185 €</b>
Digital control for micromotor DEFINITIVE	<b>376 €</b>

## Dental unit

### ABC 04-012 (2 modification)



ABC 04-012 (2modification) with lower handpieces delivery system. In set: block, saliva ejector suction system for spittle, equipment table for five equipments on pantograph arm, bottle water system, multifunctional pedal and table for instruments.	<b>1 860 €</b>
ABC 04-012 (2modification) with upper handpieces delivery system	<b>2 100 €</b>
Autonomic drain system	<b>140 €</b>
<b>Complekting operating light</b>	
Halogen dental light with two light intensity on pantograph arm	<b>223 €</b>
LED dental light DL60/DL66 with 1 LED, two light intensity on pantograph arm	<b>376 €</b>
LED dental light DL150/DL155 with 4 LED, three light intensity on pantograph arm	<b>400 €</b>
LED dental light with 6 LED, three light intensity on pantograph arm	<b>360 €</b>
<b>Variable equipment set for unit</b>	
Turbine connector	<b>133 €</b>
Turbine connector with light	<b>188 €</b>
Micromotor MV2.20 with connector	<b>160 €</b>
Air micromotor VORTICE (with turbo connection)	<b>300 €</b>
Micromotor brushless with connector DEFINITIVE SET	<b>350 €</b>
Micromotor brushless with LED, with connector DEFINITIVE LED SET	<b>435 €</b>
Technical micromotor 102 L with connector	<b>266 €</b>
Technical micromotor 108 E with connector	<b>226 €</b>
Multifunctional syringe – angular, 3 function	<b>70 €</b>
Multifunctional syringe – straight, 3 function	<b>115 €</b>
Multifunctional syringe – straight /angular, 6 function	<b>220 €</b>
Multifunctional syringe Minilight Fibra Optical – straight /angular, 7 function	<b>261 €</b>
Multifunctional syringe Minibright – straight /angular, 7 function	<b>232 €</b>
LED curing lamp	<b>480 €</b>
Ultrasonic scaler	<b>510 €</b>
Airpolishing for turbine coupling „Easyjet Pro“	<b>460 €</b>
Negotoscope	<b>130 €</b>
<b>Suction system with assistant table</b>	
Injector air suction system for dust	<b>270 €</b>
Vacuum suction system with automatic air/water separation “METASYS” (without vacuum pump)	<b>410 €</b>
<b>Additional dental equipment completing</b>	
Intra-oral camera with monitor and holder	<b>1 235 €</b>
Vacuum pump in box ABC 04-VS	<b>695 €</b>
Table for devices	<b>45 €</b>
Two metal armrests for patient chair (can be one armrest)	<b>185 €</b>
Digital control for micromotor DEFINITIVE	<b>376 €</b>

## Dental unit

### SIGER U100



SIGER U 100 chair-mounted, with lower handpieces delivery system. In set: block, patient chair ABC 04-KP with left plastic armrest, saliva ejector air suction system for spittle, equipment table for five equipments on pantograph arm, bottle water system, multifunctional pedal and table for instruments.	<b>5 130 €</b>
SIGER U 100 chair-mounted, with upper handpieces delivery system	<b>5 370 €</b>
Autonomic drain system	<b>140 €</b>
<b>Complekting operating light</b>	
Halogen dental light with two light intensity on pantograph arm	<b>223 €</b>
LED dental light DL60/DL66 with 1 LED, two light intensity on pantograph arm	<b>376 €</b>
LED dental light DL150/DL155 with 4 LED, three light intensity on pantograph arm	<b>400 €</b>
LED dental light with 6 LED, three light intensity on pantograph arm	<b>360 €</b>
<b>Variable equipment set for unit</b>	
Turbine connector	<b>133 €</b>
Turbine connector with light	<b>188 €</b>
Micromotor MV2.20 with connector	<b>160 €</b>
Air micromotor VORTICE (with turbo connection)	<b>300 €</b>
Micromotor brushless with connector DEFINITIVE SET	<b>350 €</b>
Micromotor brushless with LED, with connector DEFINITIVE LED SET	<b>435 €</b>
Technical micromotor 102 L with connector	<b>266 €</b>
Technical micromotor 108 E with connector	<b>226 €</b>
Multifunctional syringe – angular, 3 function	<b>70 €</b>
Multifunctional syringe – straight, 3 function	<b>115 €</b>
Multifunctional syringe – straight /angular, 6 function	<b>220 €</b>
Multifunctional syringe Minilight Fibra Optical – straight /angular, 7 function	<b>261 €</b>
Multifunctional syringe Minibright – straight /angular, 7 function	<b>232 €</b>
LED curing lamp	<b>480 €</b>
Ultrasonic scaler	<b>510 €</b>
Airpolishing for turbine coupling „Easyjet Pro“	<b>460 €</b>
Negotoscope	<b>130 €</b>
<b>Suction system with assistant table</b>	
Injector air suction system for dust	<b>270 €</b>
Vacuum suction system with automatic air/water separation “METASYS” (without vacuum pump)	<b>410 €</b>
<b>Additional dental equipment completing</b>	
Intra-oral camera with monitor and holder	<b>1 235 €</b>
Vacuum pump in box ABC 04-VS	<b>695 €</b>
Table for devices	<b>45 €</b>
Two metal armrests for patient chair (can be one armrest)	<b>185 €</b>
Digital control for micromotor DEFINITIVE	<b>376 €</b>

## Mobil dental unit ABC 04

MOBIL DENTAL UNIT "ABC 04" with equipment table for four equipments, bottle water system and operating pedal. **390 €**

### Variable equipment set for unit

Turbine connector	133 €
Turbine connector with light	188 €
Micromotor MV2.20 with connector	160 €
Air micromotor VORTICE (with turbo connection)	300 €
Micromotor brushless with connector DEFINITIVE SET	350 €
Micromotor brushless with LED, with connector DEFINITIVE LED SET	435 €
Technical micromotor 102 L with connector	266 €
Technical micromotor 108 E with connector	226 €
Multifunctional syringe – angular, 3 function	70 €
Multifunctional syringe – straight, 3 function	115 €
Multifunctional syringe – straight /angular, 6 function	220 €
Multifunctional syringe Minilight Fibra Optical – straight /angular, 7 function	261 €
Multifunctional syringe Minibright – straight /angular, 7 function	232 €
LED curing lamp	480 €
Ultrasonic scaler	510 €
Airpolishing for turbine coupling „Easyjet Pro“	460 €
Negotoscope	130 €



## Dental unit ABC 04-06

Dental unit "ABC 04-06" is suitable for dentist-technician, who in work not uses water. It is fix on table or other operating surface. In set: equipment table for three equipments, block and control pedal. **132 €**

### Variable equipment set for unit

Turbine connector	133 €
Micromotor MV2.20 with connector	160 €
Technical micromotor 102 L with connector	266 €
Technical micromotor 108 E with connector	226 €
Multifunctional syringe – angular, 3 function	70 €
LED curing lamp	480 €
Reducer	30 €



(with foot control)

## Micromotor ABC 04-05 K

Mobile dental micromotor with pedal for different types of micromotors handpieces **428 €**

## Dental canula ABC 04-SA

Dental canula for suction saliva. With vacuum pump, cuspidor and saliva collect system **290 €**



## Mobile dental canula ABC 04-AS

Mobile dental canula for suction dust and saliva. With vacuum pump and drain collect system **900 €**



## Mobil dental vacuum station ABC 04-VS

Vacuum pump in box, with control **695 €**

## Patient chair ABC 04 - KP

Patient chair with left plastic armrest **1 795 €**

Two metal armrests for patient chair (can be one armrest) **185 €**



OK-1

OK-2

OKR-2

## Dental chairs

Dental chair "OK-1", "OK-2" (with foot control) **245 €**

Dental chair "OKR-1", "OKR-2" (with hand control) **190 €**

Dental chair "OKR-saddler" (with hand control) **250 €**

Dental chair "OK-besides flap" (with foot control) **200 €**

Dental chair "OKR-besides flap" (with hand control) **145 €**



OKR-saddler

OKR, OK besides flap

## Dental cabinet

Dental cabinet (with trundle) **240 €**





E. Navicko enterprise "Medica"  
Architektų str. 1, LT-78366 Siauliai, Lithuania  
Tel.: +370 41 514162, +370 41 514164  
Fax +370 41 514162  
Mob. +370 676 41208  
E-mail: info@medica.lt www.medica.lt

## ***Palette of colors for dental equipment***

